

Program Your Replacement Pad/Mat to TL-2100G

Quick Start Instructions

Program Sensor Pad(s)

*Before programming: Activate each pad by applying 30 seconds of pressure.

The CordLess Monitor can work in two modes:


1. With one CordLess sensor pad
2. With two CordLess sensor pads

(See separate instructions for CordLess floor mat operation)


Place one or two CordLess bed or chair pads within reach. Each sensor pad comes with a CordLess transmitter connected and turned on.

Press and hold the program button on the back of the monitor until you hear three short beeps to enter programming mode. Release the button.

*NOTE: When there are two pads programmed to the monitor and no one is on the pads, press the reset button twice to reset the monitor.



Call Customer Service at (800) 650-3637 for assistance



Program Your Replacement Pad/Mat to TL-2100G

Quick Start Instructions

Program Sensor Pad(s)

*Before programming: Activate each pad by applying 30 seconds of pressure.

The CordLess Monitor can work in two modes:


1. With one CordLess sensor pad
2. With two CordLess sensor pads

(See separate instructions for CordLess floor mat operation)

Place one or two CordLess bed or chair pads within reach. Each sensor pad comes with a CordLess transmitter connected and turned on.

Press and hold the program button on the back of the monitor until you hear three short beeps to enter programming mode. Release the button.

*NOTE: When there are two pads programmed to the monitor and no one is on the pads, press the reset button twice to reset the monitor.



Call Customer Service at (800) 650-3637 for assistance

Transmitter Instructions

Replacement Pads/Mats

Please make sure your transmitter inside the neck of your sensor pad/mat is turned on before trying to program monitor


1. Open neck of pad and pull out transmitter.
2. Turn on the transmitter by moving the switch to the on position.
3. Once you have the switch in the on position, you should see the red light turn on and stay lit for a short time on the top side on the transmitter.
4. Once the red light turns off, it should blink every 45 seconds to show that it has power.
5. Now that you have the transmitter turned on, you can place it back inside the neck of the pad and seal closed with provided sticker. Please activate your pad following the steps below.

Pads must be activated prior to first use.
Pads are shipped in "Sleep Mode" to preserve battery life.

To activate:

Press and hold pad for 30 seconds then release. This pad is now activated and "ON". This is only necessary for the first time use.

-See monitor instructions for programming



*Once Transmitter is turned on, it is intended to stay on for the life of the pad. DO NOT use as an ON/OFF switch.

Press and hold for 30 seconds



Transmitter Instructions

Replacement Pads/Mats

Please make sure your transmitter inside the neck of your sensor pad/mat is turned on before trying to program monitor


1. Open neck of pad and pull out transmitter.
2. Turn on the transmitter by moving the switch to the on position.
3. Once you have the switch in the on position, you should see the red light turn on and stay lit for a short time on the top side on the transmitter.
4. Once the red light turns off, it should blink every 45 seconds to show that it has power.
5. Now that you have the transmitter turned on, you can place it back inside the neck of the pad and seal closed with provided sticker. Please activate your pad following the steps below.

Pads must be activated prior to first use.
Pads are shipped in "Sleep Mode" to preserve battery life.

To activate:

Press and hold pad for 30 seconds then release. This pad is now activated and "ON". This is only necessary for the first time use.

-See monitor instructions for programming



*Once Transmitter is turned on, it is intended to stay on for the life of the pad. DO NOT use as an ON/OFF switch.

Press and hold for 30 seconds

