

Always-On System Reduces Risk

Smart Caregiver's Safety-Reset Monitor is the most reliable system on the market. No on/off switch means no forgetting to turn it on and it cannot be defeated by residents. Available for both chair and bed, these monitors engage at the moment the patient sits or lays down on the pad. Caregivers know the system is engaged and operating at all times by an indicator LED. Easy for staff to use. And, because it is wireless, one alarm is outside the room where it won't disturb others. Alarm at bed can be adjusted from silent to bud. Guaranteed for one year.

- Bigger, Wider, Pad Made from Soft FDA Approved Vinyl
- Monitor Runs On AC and Battery
- An Always-On System Is Simpler and Safer
- Also Plugs Into Nurse Call System
- Microprocessor Extends Battery Life
- Also Available with Disposable 90 Day Pads

For more information, please call 800 650-3637 or log on to www.smartcaregivercorp.com