

Prevent Falls & Wandering

Easy-To-Install Call Button & Light System



**Arrives pre-programmed
and ready to use!**



2007-SYSR1

Easy Installation

The 2007-SYSR1 Emergency Call Light System includes a wireless nurse call button with pull cord, reset button, and call alarm with light and is ideal for bathroom and hallway notification systems. All components are easily mounted on most surfaces with no hard wiring required.

Key Features:

- Waterproof!
- Audible & Visible Alert
- No Hard Wiring
- Can be programmed to work with the TL-2015 Large Facility Central Monitor
- Self-Adhesive Wall Mounting
- Ideal for Bathroom Use



Call Us Today! (800) 650-3637 | www.smartcaregiver.com