


# Prevent Falls & Wandering

## Anti-Wandering Door Systems

Available in single or double door systems.

For facilities of any size!



### All-in-One

This "All in One" Anti-Wandering Door Bar System functions with just a door alarm and a resident wristband. Simply mount the door monitor alarm by any exit or doorway in your care facility, plug it in, and place wristband on resident. When a resident wearing a wristband attempts to wander too close or through the doorway, the door monitor alarm sounds audibly and visually.

### Simple Installation

1. Mount door bar & activate transmitter
2. Test
3. Put wristband on resident

**Affordable – No Wiring – No Electricians**

MFR#	PRODUCT DESCRIPTION
TL-3005SYSR2	Single Door Bar Monitor System (1 year warranty)
TL-3004SYSR2	Double Door Bar Monitor System (1 year warranty)
TL-2012SR1	Resident Wristband Transmitter (6 month warranty)
TL-2014W	Wall Mount Reset Button (1 year warranty)
TL-2015R2	Large Facility Central Monitoring Unit (1 year warranty)
DT-01	Door Tester & Wristband Control Unit (1 year warranty)
TL-2025	Magnetic Door Lock (1 year warranty)

Also works with:



**Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)**