


Lights Turn On When They Get Up

No getting up in the dark!



- Requires Paired Sensor Device

Eliminate Nighttime Obstacles

The Smart Light Outlet will turn on a light when triggered by the Cordless Bed Sensor Pad or floor mat. The Smart Light Outlet allows you to plug in any standard lamp or light source to create a safer environment by providing illumination when someone gets out of bed. The light will remain on until pressure is placed back on the pad or the outlet has been reset by the caregiver.

Use With Any Smart CordLess® Pads or Mats

The Smart Outlet must be programmed to a CordLess Sensor Device which is then paired to a CordLess monitor, which can be placed remotely to alert caregivers at a designated location. The Smart Light Outlet will turn on a light when pressure is removed from a bed or chair sensor pad or when pressure is placed on a floor mat or landing mat.



SMS-11 Smart Light Socket

Key Features:

- Immediately illuminates room upon activation!
- Cord Free Operation
- Turns lights on before residents gets up
- Works with sensor pads and floor mats



Fall Prevention by Smart Caregiver
(800) 650-3637 | www.smartcaregiver.com