

# Prevent Falls & Wandering

## Corded Pads & Mats For Bed & Chair Exit Alarms



### Works with:



Wireless Fall Monitor  
TL-2016



Basic Fall Monitor  
TL-2100B



Economy Fall Monitor  
TL-2100E



Safety-Reset Fall Monitor  
TL-2100S



Recordable Voice Monitor  
TL-3100V

### Simple early warning system to help prevent falls.

Bed and chair sensor pads alert caregivers whenever pressure is removed from the sensor pad (i.e. when a resident gets up from a bed or chair & pressure is removed from the pad).

Floor Mats & Landing Mats alert caregivers when pressure is applied to the sensor mat (i.e. when a resident steps onto a floor mat and pressure is placed on the mat).

Floor mats are ideal for bedside and doorway mobility monitoring and bed and chair pads are ideal early warning alerts of resident mobility as they exit beds and chairs.

MFR#	SMART PRODUCT DESCRIPTION
PPB-WI	Bed Pad 20" x 30" (1 year pad life)
PPB-RI	Bed Pad 10" x 30" (1 year pad life)
PPB-90R, PPB-45R	Bed Pad 10" x 30" (90 day, 45 day pad life)
PPC-WI	Chair Pad 10" x 15" (1 year pad life)
PPC-90, PPC-45	Chair Pad 10" x 15" (90 day, 45 day pad life)
FM-03	Floor Mat 24" x 36" Black (1 year warranty)
FM-05	Floor Mat 24" x 36" Grey (1 year warranty)
FM-07	Floor Mat 24" x 48" Grey (1 year warranty)
LM-01	Corded Landing Mat 24" x 48", Grey (1 year warranty)



Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)