


# Prevent Falls & Wandering

## Wireless Central Monitoring



## Complete Resident Monitoring & Paging System

Wireless Deluxe Central Monitoring Unit – Designed to wirelessly connect with up to 999 devices. Monitor alerts caregiver when a resident is in need of assistance and displays which device triggered the alarm. The system can identify devices by room number.

The Deluxe CMU has a 1 year warranty and functions with a variety of components and accessories including Anti-Wandering Door Systems & LCD Pagers.

Easy setup and works with a large selection of wireless devices without any hard-wiring.

### Key Features:

- Monitor as many residents as needed (up to 999)
- Alerts indicate which device triggered the alarm
- Sends alert to caregiver pager
- Receives signal from any number of wireless devices in the facility (i.e. Anti-Wandering Door Bars, Wireless Fall Monitors, Cordless Bed Pads, Floor mats, Emergency Call Buttons, Bathroom Call Alarms, etc.)



Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)