

Prevent Falls & Wandering

Motion Activated Caregiver Paging System



Easy-To-Use Caregiver Alert!

No Noise in Room!

The Smart Caregiver TL-5102MP is a stand-alone motion paging system that allows you to eliminate in-room alarm noise. The Motion Sensor is placed in a doorway or by the bedside. When motion is detected, a wireless signal is sent to the caregiver pager up to 300 feet away.

Caregiver Alerted by Pager

The Caregiver Pager is pocket-sized and light-weight; giving caregivers the freedom they need.

Place it where you need it

The sensor can be placed next to the bed, by a doorway or anywhere to detect movement.

TL-5102MP

Key Features:

- Quiet - No alarming monitors
- Easy to set up - insert batteries and it's ready to use
- Pocket sized pager - Easy to carry
- Motion Sensor can be placed by bed or doorway to monitor residents activities
- Batteries included (2-AA, 3-AA)
- Optional AC adapter available (AC-04)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com