

Prevent Falls & Wandering

NEW
Available Now!

Facility Grade Large Central Monitor

Easy setup and works with large selection of wireless devices without any hard-wiring!



DepartAlert® Door Bar System
TL-4005SYS

Large Facility Central Monitor
TL-4015



Caregiver Pager
TL-4015P



Quiet Wireless Monitor
TL-4016R



Soft Touch Button
NC-03

Reset Button
4017-RB



Call Button with Emergency Pull Cord
4017-CB



Emergency Call Alert
4017-ECA



Corded Bed Pad



Corded Chair Pad



Corded Floor Mat

TL-4015 Key Features:

- Alert indicates which device triggered the alarm
- Sends alert to caregiver pager
- Receives signal from up to 99 wireless devices in the facility
- Batteries not included (4 x C)
- AC Adapter included (AC-04)
- Range: up to 300 feet

Don't Let Them Walk Alone

Designed to wirelessly connect with up to 99 devices. Alerts caregiver when resident is in need of assistance and displays which device triggered the alarm. The Smart Station can identify alerting devices with a designated number. The Smart Station functions with a variety of components and accessories including, Anti-Wandering Door Systems, Emergency Call Button & Light Systems, LCD Pagers and Quiet Fall Prevention Alert to Pager Systems.

Call Customer Service at (800) 650-3637 or visit www.smartcaregiver.com

