


Prevent Falls & Wandering

Smart Wireless Caregiver Personal Pager System



Monitor multiple residents from one central location

Eliminate In-Room Alarm Noise

The Smart Caregiver 433-CMU Central Monitoring Unit allows you to monitor multiple wireless components from one central location, such as a Nurses' Station. This eliminates in-room alarm noise & significantly reduces the tripping hazards associated with corded units.

Know Who Needs Help

The Central Monitor alerts audibly and the display visually indicates which resident needs assistance. At a glance, this provides caregivers with the information they need to protect their residents from falls.

Local Pager Notification

The CMU communicates with any number of 433-Pager Units, and upon activation, the pager audibly alerts and LCD screen displays which component activated the alert. This enables caregivers to attend to their duties yet be alerted in real time that a resident needs assistance. 433-PRB has a reset button that resets the 433-CMU-40/60 when alarming.

Key Features:

- Alarms audibly & visually at a central location and portable caregiver pager(s)
- Compatible with all Smart wireless components
- 40 & 60 channel capacity models available
- Range up to 300 feet
- Systems may be pre-programmed prior to delivery or easily re-programmed using the provided instructions
- AC adapter included (AC-04)
- Batteries not included (4-C)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com