

# Prevent Falls & Wandering

## Reduce Falls at Night Turn the lights on when they get up!



Smart Outlet

### Use With Smart CordLess Pads or Floor Mats

The Smart Outlet will turn on a light when pressure is removed from a bed or chair sensor pad or when pressure is placed on a floor mat or landing mat.

### Lights Turn On To Reduce Tripping Hazards

The Smart Outlet will turn on a light when triggered by the Cordless Bed Sensor Pad / Floor Mat. The Smart Outlet allows you to plug in any standard lamp or light source to create a safer environment by providing light when someone gets out of bed.

### Works With Remote Receivers

The receivers can be placed remotely to alert caregivers at a designated location.

### Key Features:

- Immediately lights up room upon activation!
- Cord Free Operation
- Turns lights on when resident gets up
- Works with CordLess sensor pads and floor mats

**\*Must be programmed with a CordLess Sensor Pad / Floor Mat.**



Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)