Prevent Falls & Wandering

Easy-To-Use **Economy Bed / Chair Exit Alarm**





Affordable Easy-to-Use Fall Monitor

This Smart Caregiver Economy Fall Monitor Alarm functions with weight-sensing bed and chair pressure pads, floor mats, or early warning seat belts. Caregiver is notified by an audible and visual alert.

Monitor can be placed at the bedside or wheelchair and then connected to a corded sensor pad, floor mat or seat belt. When the resident is about to get up from a bed or a chair, the monitor will sound alerting the caregiver.

On/Off Switch

The fall monitor features an "On/Off" switch that prevents the alarm from sounding when the caregiver moves the resident from bed or chair, or when changing clothing or linens.

TL-2100E

Key Features:

- Full-Range Volume Control
- Nurse Call Capability enables you to use with existing hard-wired nurse call systems
- Low Battery & Status Indicator Lights
- Batteries not included (3-AA)

