

Prevent Falls & Wandering

Motion Sensing Exit Alarm



Stand-Alone Monitoring Solution

The Smart Caregiver TL-2700 allows you to maintain effective awareness of resident movement without the use of pressure pads, floor mats or pull-string type devices.

Use It Where You Need It

Motion sensor can be placed at the bedside to alert when someone attempts to get out of bed, by a chair to alert when someone attempts to get up from a chair, or by a doorway to alert when someone attempts to exit the area. Bracket may be wall-mounted with provided hardware and optional mounting bracket accessories are available.

Wide Area of Detection

This motion sensing exit alarm has a wide field of detection that extends outwards approximately 30 feet from the alarm. Any movement that interrupts this field will trigger the alarm.

TL-2700



Key Features:

- Chime Alert
- Great for bedsides, hallways or doorways
- Stand-alone motion sensor alarm
- No cords or wires required
- Battery included (1-9V)
- Manual Off/Reset



Call Us Today! (800) 650-3637 | www.smartcaregiver.com