

Don't Let Them Walk Away!


New & Affordable!

Prevent Wandering With Anti-Wandering Door Systems

For Homes & Facilities of any size!



Wireless Signal



Affordable – No Wiring – No Electricians

MFR#	PRODUCT DESCRIPTION
TL-4005SYS	Single Door Bar Monitor System (1 year warranty)
TL-4012	Resident Wristband Transmitter (6 month warranty)
TL-4014W	Wall Mount Reset Button (1 year warranty)

Optional Additional Components:



TL-2025
Magnetic
Door Lock



TL-4015
Central Monitor

All-in-One

This "All in One" Anti-Wandering Door Bar System functions with just a door alert and a resident wristband. Simply mount the door bar by any exit or doorway in your care facility, plug it in, and place wristband on resident. When a resident wearing a wristband gets within 8 to 14 feet of the door bar, the door bar will alert audibly and visually.

Simple Installation

1. Mount Door Bar
2. Activate Wristband Transmitter
3. Test



Authorized Distributor
SMART
Quieter Fall Prevention

Call Us Today! (800) 650-3637 | www.smartcaregiver.com