

CMS  
Compliant Fall  
Prevention

# Prevent Falls & Wandering

## Smart Wireless Caregiver Personal Pager™ Removes noise from room



### Reduce Alarm Noise in Room!

The monitor can be set to silent to reduce in-room alarm noise. When a resident attempts to get up, a silent wireless signal is sent to a caregiver pager or central nurse's station.

### Flexible & Versatile

The wireless fall monitor is designed to work with Corded Bed & Chair Sensor Pads, Floor Mats, Seat Belts, Pagers, Emergency Call Lights, Nurse Call Cords and Central Monitoring Units. This monitor allows caregivers the flexibility to monitor residents based on their unique needs.

### Increase Caregiver Freedom!

Caregiver can attend to other duties and receive an alert with use of a pager when someone is getting up from a bed or chair. Can also be used with a floor mat to notify if they are approaching a door.

Quiet  
**WIRELESS™**

TL-4016



CareGiver  
**PAGER™**

TL-4016P

### Key Features:

- Monitor sends wireless signal to pager up to 300 feet away
- Call Button on front of monitor
- Nurse Call Capability - enables you to use with existing hard wired nurse call system
- Batteries not included (3-AA)
- Optional AC adapter available (AC-05)
- Works with Large Facility  
TL-4015 CMU

Anti-Wandering Solutions  
**SMART®**  
Quieter Fall Prevention™

Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)