

Prevent Falls & Wandering

Know When They Get Up! Low Cost Fall Prevention Monitor



TL-2100B

Great Value! All The Functions

The Smart Caregiver TL-2100B Monitor is a Full-Function, Low Cost fall monitor designed for use with all Smart Caregiver Corded Pads, Floor Mats and Seat Belts. This is a high quality product at a great, everyday value.

Auto-Reset Function

Monitor can be placed at the bedside or wheelchair, connected to a corded sensor pad. The sensor pad is placed under the resident. When the resident gets up and pressure is removed from the pad, the monitor will sound alerting the caregiver. When pressure is re-applied to the pad, the monitor will silence and reset for use.

Monitor Can Be Set To Silent

Monitor can be set to silent only when using an existing nurse call system and additional nurse call cord.

Key Features:

- Full-Range Volume Control
- Nurse Call Capability - enables you to use with existing hard-wired nurse call systems
- Low Battery & Status Indicator Lights
- Battery not included (1-9V)
- Optional AC adapter available (AC-02)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com